

# R Packages for Noncompartmental Analysis

*Kyun-Seop Bae*

*Sungpil Han*

*2017-01-26*



# Contents

<b>Preface</b>	<b>5</b>
<b>1 Introduction</b>	<b>7</b>
<b>2 Comparison of Noncompartmental analysis software</b>	<b>9</b>
2.1 Certara Phoenix WinNonLin . . . . .	9
2.2 R package: PKNCA . . . . .	9
2.3 R package: ncar . . . . .	10
2.4 R package: ncappc . . . . .	10
2.5 R package: PK . . . . .	10
<b>3 Details</b>	<b>11</b>
3.1 ncar . . . . .	11
3.2 PKNCA: Automation of Noncompartmental Analysis in R . . . . .	14



# Preface

In this book, we will cover R Packages for Noncompartmental Analysis.



# Chapter 1

## Introduction

There are growing number of open source noncompartmental analysis R packages and we will compare, review and discuss them.





## Chapter 2

# Comparison of Noncompartmental analysis software

We will perform NCA with a dataset **Theoph** and compare  $C_{\max}$  and  $AUC_{\inf}$  calculated by each software or package.

### 2.1 Certara Phoenix WinNonLin

<https://www.certara.com/software/pkpd-modeling-and-simulation/phoenix-winnonlin/>

#### 2.1.0.1 Pros

- Validated for several years
- Industry standard
- Versatile unit setting
- Easy using by GUI
- Generating plots supported

#### 2.1.0.2 Cons

- Expansive (~several thousand dollars)
- Not suitable for reproducible research
- CDISC SDTM not compatible (input and output)

### 2.2 R package: PKNCA

<https://github.com/billdenney/pknca>

#### 2.2.0.1 Pros

- Open source and free of charge
- CDISC SDTM semi compatible (output)
- Calculate partial(interval) AUC with ‘linear’ or ‘log’ interpolation method but in a cumbersome way

### 2.2.0.2 Cons

- CDISC SDTM not compatible (input)
- More tests required
- Experience with R language required
- Generating plots not supported for now (To be supported soon)

## 2.3 R package: ncar

Extension of NonCompart for CDISC input

<https://github.com/asancpt/ncar>

### 2.3.0.1 Pros

- Open source and free of charge
- CDISC SDTM fully compatible (input - DM, EX, PC)
- Output of CDISC SDTM PP domain terms
- Calculate partial(interval) AUC with ‘linear’ or ‘log’ interpolation method

### 2.3.0.2 Cons

- More tests required
- Experience with R language required
- Generating plots not supported for now (To be supported soon)

## 2.4 R package: ncappc

NCA Calculation and Population PK Model Diagnosis

<https://cran.r-project.org/web/packages/ncappc/index.html>  
27000291

<https://www.ncbi.nlm.nih.gov/pubmed/>

```
#install.packages("ncappc")
#library(ncappc) # Broken Link
```

## 2.5 R package: PK

Basic Non-Compartmental Pharmacokinetics

<https://cran.r-project.org/web/packages/PK/index.html>

```
#install.packages("PK")
library(PK)
```

```
## ***** PK Version 1.3-3*****
```

```
## Type PKNews() to see new features/changes/bug fixes.
```

# Chapter 3

## Details

### 3.1 ncar

```
library(ncar)
```

```
NCA(Theoph, "Subject", "Time", "conc", Dose=320)
```

##	Subject	C <sub>MAX</sub>	C <sub>MAXD</sub>	T <sub>MAX</sub>	T <sub>LAG</sub>	CL <sub>ST</sub>	CL <sub>STP</sub>	TL <sub>ST</sub>	LAM <sub>ZHL</sub>
## 1	1	10.50	0.03281250	1.12	0	3.28	3.2801465	24.37	14.304378
## 2	2	8.33	0.02603125	1.92	0	0.90	0.8886398	24.30	6.659342
## 3	3	8.20	0.02562500	1.02	0	1.05	1.0550967	24.17	6.766087
## 4	4	8.60	0.02687500	1.07	0	1.15	1.1564216	24.65	6.981247
## 5	5	11.40	0.03562500	1.00	0	1.57	1.5556951	24.35	8.002264
## 6	6	6.44	0.02012500	1.15	0	0.92	0.9412712	23.85	7.894998
## 7	7	7.09	0.02215625	3.48	0	1.15	1.1607192	24.22	7.846668
## 8	8	7.56	0.02362500	2.02	0	1.25	1.2285268	24.12	8.510038
## 9	9	9.03	0.02821875	0.63	0	1.12	1.1164831	24.43	8.405999
## 10	10	10.21	0.03190625	3.55	0	2.42	2.4136923	23.70	9.246916
## 11	11	8.00	0.02500000	0.98	0	0.86	0.8598066	24.08	7.261237
## 12	12	9.75	0.03046875	3.52	0	1.17	1.1755390	24.15	6.286508
##		LAM <sub>Z</sub>	LAM <sub>ZLL</sub>	LAM <sub>ZUL</sub>	LAM <sub>ZNPT</sub>	COR <sub>RY</sub>	R <sup>2</sup>	R <sup>2</sup> <sub>ADJ</sub>	
## 1	0.04845700	9.05	24.37		3	-0.9999999	0.9999997	0.9999995	
## 2	0.10408644	7.03	24.30		4	-0.9985967	0.9971954	0.9957931	
## 3	0.10244431	9.00	24.17		3	-0.9996624	0.9993250	0.9986499	
## 4	0.09928702	9.02	24.65		3	-0.9994619	0.9989241	0.9978483	
## 5	0.08661888	7.02	24.35		4	-0.9993234	0.9986472	0.9979708	
## 6	0.08779574	2.03	23.85		7	-0.9991203	0.9982413	0.9978896	
## 7	0.08833650	6.98	24.22		4	-0.9993349	0.9986702	0.9980053	
## 8	0.08145054	3.53	24.12		6	-0.9954961	0.9910124	0.9887655	
## 9	0.08245863	8.80	24.43		3	-0.9997218	0.9994437	0.9988873	
## 10	0.07495982	9.38	23.70		3	-0.9997543	0.9995087	0.9990174	
## 11	0.09545856	9.03	24.08		3	-0.9999991	0.9999983	0.9999965	
## 12	0.11025949	9.03	24.15		3	-0.9996984	0.9993968	0.9987936	
##		AUC <sub>LST</sub>	AUC <sub>ALL</sub>	AUC <sub>IFO</sub>	AUC <sub>IFOD</sub>	AUC <sub>PEO</sub>	AUC <sub>IFP</sub>	AUC <sub>IFPD</sub>	
## 1	148.92305	148.92305	216.61193	0.6769123	31.248917	216.61496	0.6769217		
## 2	91.52680	91.52680	100.17346	0.3130421	8.631687	100.06432	0.3127010		
## 3	99.28650	99.28650	109.53597	0.3422999	9.357173	109.58572	0.3424554		
## 4	106.79630	106.79630	118.37888	0.3699340	9.784331	118.44356	0.3701361		

```
## 5 121.29440 121.29440 139.41978 0.4356868 13.000579 139.25463 0.4351707
## 6 73.77555 73.77555 84.25442 0.2632951 12.437174 84.49670 0.2640522
## 7 90.75340 90.75340 103.77180 0.3242869 12.545221 103.89315 0.3246661
## 8 88.55995 88.55995 103.90669 0.3247084 14.769730 103.64305 0.3238845
## 9 86.32615 86.32615 99.90872 0.3122147 13.594978 99.86607 0.3120815
## 10 138.36810 138.36810 170.65206 0.5332877 18.918002 170.56791 0.5330247
## 11 80.09360 80.09360 89.10274 0.2784461 10.110962 89.10072 0.2784397
## 12 119.97750 119.97750 130.58883 0.4080901 8.125757 130.63907 0.4082471
##      AUCPEP      AUMCLST      AUMCIFO      AUMCPEO      AUMCIFP      AUMCPEP      MRTEVLST
## 1 31.249876 1459.0711 4505.5348 67.61603 4505.6709 67.61701 9.797483
## 2 8.532030 706.5866 999.7723 29.32525 996.0716 29.06267 7.719996
## 3 9.398325 803.1859 1150.9648 30.21629 1152.6529 30.31850 8.089578
## 4 9.833594 901.0842 1303.2524 30.85881 1305.4981 30.97775 8.437410
## 5 12.897403 1017.1143 1667.7216 39.01174 1661.7937 38.79419 8.385501
## 6 12.688246 609.1524 978.4285 37.74176 986.9665 38.28034 8.256833
## 7 12.647366 782.4199 1245.0984 37.16000 1249.4111 37.37691 8.621383
## 8 14.552931 739.5346 1298.1158 43.03015 1288.5201 42.60589 8.350666
## 9 13.558076 705.2296 1201.7715 41.31750 1200.2124 41.24126 8.169363
## 10 18.878001 1278.1800 2473.9934 48.33535 2470.8765 48.27018 9.237534
## 11 10.108918 617.2422 928.5600 33.52694 928.4900 33.52193 7.706511
## 12 8.161087 977.8807 1330.3840 26.49636 1332.0528 26.58844 8.150534
##      MRTEVIFO      MRTEVIFP      VZFO      VZFP      CLFO      CLFP
## 1 20.800031 20.800368 30.48675 30.48632 1.477296 1.477276
## 2 9.980411 9.954313 30.69044 30.72392 3.194459 3.197943
## 3 10.507642 10.518276 28.51710 28.50415 2.921415 2.920088
## 4 11.009163 11.022112 27.22596 27.21110 2.703185 2.701709
## 5 11.961873 11.933490 26.49799 26.52942 2.295227 2.297949
## 6 11.612785 11.680533 43.25973 43.13569 3.798020 3.787130
## 7 11.998427 12.025924 34.90844 34.86767 3.083689 3.080088
## 8 12.493092 12.432287 37.81051 37.90669 3.079686 3.087520
## 9 12.028695 12.018220 38.84279 38.85938 3.202924 3.204292
## 10 14.497296 14.486174 25.01554 25.02788 1.875160 1.876086
## 11 10.421227 10.420679 37.62219 37.62304 3.591360 3.591441
## 12 10.187579 10.196436 22.22429 22.21575 2.450439 2.449497
```

```
iAUC = data.frame(Name=c("AUC[0-12h]", "AUC[0-24h]"), Start=c(0,0), End=c(12,24)) ; iAUC
```

```
##      Name Start End
## 1 AUC[0-12h]      0 12
## 2 AUC[0-24h]      0 24
```

```
NCA(Theoph, "Subject", "Time", "conc", Dose=320, iAUC=iAUC)
```

```
##      Subject CMAX      CMAXD TMAX TLAG CLST      CLSTP TLST      LAMZHL
## 1      1 10.50 0.03281250 1.12      0 3.28 3.2801465 24.37 14.304378
## 2      2 8.33 0.02603125 1.92      0 0.90 0.8886398 24.30 6.659342
## 3      3 8.20 0.02562500 1.02      0 1.05 1.0550967 24.17 6.766087
## 4      4 8.60 0.02687500 1.07      0 1.15 1.1564216 24.65 6.981247
## 5      5 11.40 0.03562500 1.00      0 1.57 1.5556951 24.35 8.002264
## 6      6 6.44 0.02012500 1.15      0 0.92 0.9412712 23.85 7.894998
## 7      7 7.09 0.02215625 3.48      0 1.15 1.1607192 24.22 7.846668
## 8      8 7.56 0.02362500 2.02      0 1.25 1.2285268 24.12 8.510038
## 9      9 9.03 0.02821875 0.63      0 1.12 1.1164831 24.43 8.405999
## 10     10 10.21 0.03190625 3.55      0 2.42 2.4136923 23.70 9.246916
## 11     11 8.00 0.02500000 0.98      0 0.86 0.8598066 24.08 7.261237
```

```

## 12      12  9.75 0.03046875 3.52    0 1.17 1.1755390 24.15  6.286508
##          LAMZ LAMZLL LAMZUL LAMZNPT      CORRXY      R2      R2ADJ
## 1  0.04845700   9.05  24.37      3 -0.9999999 0.9999997 0.9999995
## 2  0.10408644   7.03  24.30      4 -0.9985967 0.9971954 0.9957931
## 3  0.10244431   9.00  24.17      3 -0.9996624 0.9993250 0.9986499
## 4  0.09928702   9.02  24.65      3 -0.9994619 0.9989241 0.9978483
## 5  0.08661888   7.02  24.35      4 -0.9993234 0.9986472 0.9979708
## 6  0.08779574   2.03  23.85      7 -0.9991203 0.9982413 0.9978896
## 7  0.08833650   6.98  24.22      4 -0.9993349 0.9986702 0.9980053
## 8  0.08145054   3.53  24.12      6 -0.9954961 0.9910124 0.9887655
## 9  0.08245863   8.80  24.43      3 -0.9997218 0.9994437 0.9988873
## 10 0.07495982   9.38  23.70      3 -0.9997543 0.9995087 0.9990174
## 11 0.09545856   9.03  24.08      3 -0.9999991 0.9999983 0.9999965
## 12 0.11025949   9.03  24.15      3 -0.9996984 0.9993968 0.9987936
##          AUCCLST   AUCALL   AUCIFO   AUCIFOD   AUCPEO   AUCIFP   AUCIFPD
## 1  148.92305 148.92305 216.61193 0.6769123 31.248917 216.61496 0.6769217
## 2   91.52680  91.52680 100.17346 0.3130421  8.631687 100.06432 0.3127010
## 3   99.28650  99.28650 109.53597 0.3422999  9.357173 109.58572 0.3424554
## 4  106.79630 106.79630 118.37888 0.3699340  9.784331 118.44356 0.3701361
## 5  121.29440 121.29440 139.41978 0.4356868 13.000579 139.25463 0.4351707
## 6   73.77555  73.77555  84.25442 0.2632951 12.437174  84.49670 0.2640522
## 7   90.75340  90.75340 103.77180 0.3242869 12.545221 103.89315 0.3246661
## 8   88.55995  88.55995 103.90669 0.3247084 14.769730 103.64305 0.3238845
## 9   86.32615  86.32615  99.90872 0.3122147 13.594978  99.86607 0.3120815
## 10 138.36810 138.36810 170.65206 0.5332877 18.918002 170.56791 0.5330247
## 11  80.09360  80.09360  89.10274 0.2784461 10.110962  89.10072 0.2784397
## 12 119.97750 119.97750 130.58883 0.4080901  8.125757 130.63907 0.4082471
##          AUCPEP   AUMCLST   AUMCIFO   AUMCPEO   AUMCIFP   AUMCPEP MRTEVLST
## 1  31.249876 1459.0711 4505.5348 67.61603 4505.6709 67.61701 9.797483
## 2   8.532030  706.5866  999.7723 29.32525  996.0716 29.06267 7.719996
## 3   9.398325  803.1859 1150.9648 30.21629 1152.6529 30.31850 8.089578
## 4   9.833594  901.0842 1303.2524 30.85881 1305.4981 30.97775 8.437410
## 5  12.897403 1017.1143 1667.7216 39.01174 1661.7937 38.79419 8.385501
## 6  12.688246  609.1524  978.4285 37.74176  986.9665 38.28034 8.256833
## 7  12.647366  782.4199 1245.0984 37.16000 1249.4111 37.37691 8.621383
## 8  14.552931  739.5346 1298.1158 43.03015 1288.5201 42.60589 8.350666
## 9  13.558076  705.2296 1201.7715 41.31750 1200.2124 41.24126 8.169363
## 10 18.878001 1278.1800 2473.9934 48.33535 2470.8765 48.27018 9.237534
## 11 10.108918  617.2422  928.5600 33.52694  928.4900 33.52193 7.706511
## 12  8.161087  977.8807 1330.3840 26.49636 1332.0528 26.58844 8.150534
##          MRTEVIFO MRTEVIFP   VZFO   VZFP   CLFO   CLFP AUC[0-12h]
## 1  20.800031 20.800368 30.48675 30.48632 1.477296 1.477276  91.73552
## 2   9.980411  9.954313 30.69044 30.72392 3.194459 3.197943  67.48030
## 3  10.507642 10.518276 28.51710 28.50415 2.921415 2.920088  70.17971
## 4  11.009163 11.022112 27.22596 27.21110 2.703185 2.701709  73.05115
## 5  11.961873 11.933490 26.49799 26.52942 2.295227 2.297949  84.61490
## 6  11.612785 11.680533 43.25973 43.13569 3.798020 3.787130  51.75887
## 7  11.998427 12.025924 34.90844 34.86767 3.083689 3.080088  62.09875
## 8  12.493092 12.432287 37.81051 37.90669 3.079686 3.087520  62.71486
## 9  12.028695 12.018220 38.84279 38.85938 3.202924 3.204292  60.12123
## 10 14.497296 14.486174 25.01554 25.02788 1.875160 1.876086  90.81742
## 11 10.421227 10.420679 37.62219 37.62304 3.591360 3.591441  58.53963
## 12 10.187579 10.196436 22.22429 22.21575 2.450439 2.449497  85.02136
##          AUC[0-24h]

```

```
## 1 147.69459
## 2 91.24908
## 3 99.10481
## 4 105.99811
## 5 120.73101
## 6 73.91422
## 7 90.49567
## 8 88.40890
## 9 85.82985
## 10 139.08507
## 11 80.02431
## 12 119.79884
```

## 3.2 PKNCA: Automation of Noncompartmental Analysis in R

### 3.2.1 ISoP Pharmacometrics Study Group Presentation

- 강의 동영상 <https://www.youtube.com/watch?v=WCmFrheYtcc>
- 프로젝트 <https://github.com/bilddenney/pknca>
- Package <https://cran.r-project.org/web/packages/PKNCA/>
  - 예제 R Markdown 파일 : <https://github.com/bilddenney/pknca/tree/master/vignettes>
- PPT 파일
- PKNCA 패키지란 무엇인가? \* Pharmacokinetic(PK) data를 위한 모든 noncompartmental analysis (NCA) 계산이 가능한 R용 패키지

```
library(devtools)
install_github("bilddenney/pknca")
```

### 3.2.2 오픈소스 NCA - 지금이 적기이다.

- Data standards 가 점점 많아짐
- CDISC/SDTM가 FDA requirement
- CDISC ADaM working group is standardizing NCA data set (ADNCA) |
  - CDISC SDTM pharmacokinetic concentration (PC) and pharmacokinetic parameter (PP) domains have been standardized
- 우리도 R로 NCA?

### 3.2.3 할수 있는 것

- Organizes concentration/time and dose/time data
- Predicts what you most likely need from NCA parameters from the concentration and dosing data.
- Allows user control of all NCA parameter and summary calculations
- Calculates all (standard) NCA parameters (Targeting the SDTM PK 파라미터)
- Summarizes the parameters

### 3.2.4 한계

- 그래픽 못그림
- 파라미터의 statistics 못구함 (곧 기능 추가할듯)

### 3.2.5 PKNCA 현재는 0.7

- NCA 파라미터 계산가능 (Cmax, Tmax, AUClast, AUCinf, AUMC, half-life, ...)
- NCA-related calculations (Superposition, Concentration interpolation/extrapolation (with AUC methods), Time to steady-state)
- SDTM PP-READY OUTPUT 가능
- 인풋에서 아웃풋까지 TRACK가능하다.
- 800개 넘는 테스트 케이스가 있음.

### 3.2.6 PKNCA 곧 1.0이 나올것이다.

- Improved prediction of desired parameters (정확도 accuracy, number 등)

### 3.2.7 참고사항

- Github에서 모두 다운로드 가능
- CRAN에 package올라왔다. (0.7) <https://cran.r-project.org/web/packages/PKNCA/>  
– wdenney@humanpredictions.com 으로 메일 보내라
- 모든게 오픈이기 때문에 Github에서 기여 환영

### 3.2.8 RStudio를 사용한 Hands-on 실습

#### 3.2.8.1 Example-theophylline.Rmd

- Theophylline 농도를 가지고 PK Parameter 구하는 법
- <https://raw.githubusercontent.com/bilddenney/pknca/master/vignettes/Example-theophylline.Rmd>
- 이 파일을 RStudio에서 실행해본다.
- 이후 article에서 분석할 것입니다.

#### 3.2.8.2 Superposition.Rmd

- <https://raw.githubusercontent.com/bilddenney/pknca/master/vignettes/Superposition.Rmd>
- 이 파일을 RStudio에서 실행해본다.

### 3.2.9 Closing

- PKNCA.options() 모든 옵션을 볼 수 있다.

### 3.2.10 결론

- 써보고 feedback주고 contribute해라.

```
## It is always a good idea to look at the data
knitr::kable(head(datasets::Theoph))
```

Table 3.1: Example dosing data extracted from theophylline data set

	Dose	Time	Subject
1	4.02	0	1
12	4.40	0	2
23	4.53	0	3
34	4.40	0	4
45	5.86	0	5
56	4.00	0	6
67	4.95	0	7
78	4.53	0	8
89	3.10	0	9
100	5.50	0	10
111	4.92	0	11
122	5.30	0	12

Subject	Wt	Dose	Time	conc
1	79.6	4.02	0.00	0.74
1	79.6	4.02	0.25	2.84
1	79.6	4.02	0.57	6.57
1	79.6	4.02	1.12	10.50
1	79.6	4.02	2.02	9.66
1	79.6	4.02	3.82	8.58

```
## By default it is groupedData; convert it to a data frame for use
my.conc <- PKNCAdconc(as.data.frame(datasets::Theoph), conc~Time|Subject)
```

```
## Dosing data needs to only have one row per dose, so subset for
## that first.
```

```
d.dose <- unique(datasets::Theoph[datasets::Theoph$Time == 0,
                                c("Dose", "Time", "Subject")])
```

```
knitr::kable(d.dose,
              caption="Example dosing data extracted from theophylline data set")
```

```
my.dose <- PKNCAdose(d.dose, Dose~Time|Subject)
```

```
my.data.automatic <- PKNCAdata(my.conc, my.dose)
knitr::kable(PKNCAdoptions("single.dose.aucs"))
```

start	end	auclast	aucall	aumclast	aumcall	cmax	cmin	tmax	tlast	tfirst	clast.obs	f
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS

```
knitr::kable(my.data.automatic$intervals)
```



start	end	auclast	aucall	aumclast	aumcall	cmax	cmin	tmax	tlast	tfirst	clast.obs	f
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS
0	24	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALS
0	Inf	FALSE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALS

```
my.intervals <- data.frame(start=0,
                           end=Inf,
                           cmax=TRUE,
                           tmax=TRUE,
                           aucinf.obs=TRUE,
                           auclast=TRUE)
my.data.manual <- PKNCAdata(my.conc, my.dose,
                           intervals=my.intervals)
knitr::kable(my.data.manual$intervals)
```

start	end	auclast	aucall	aumclast	aumcall	cmax	cmin	tmax	tlast	tfirst	clast.obs	f
0	Inf	TRUE	FALSE	FALSE	FALSE	TRUE	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE

```
my.results.automatic <- pk.nca(my.data.automatic)
knitr::kable(head(my.results.automatic$result))
```

start	end	Subject	PPTESTCD	PPORRES
0	24	1	auclast	92.3654416
0	Inf	1	cmax	10.5000000
0	Inf	1	tmax	1.1200000
0	Inf	1	tlast	24.3700000
0	Inf	1	lambda.z	0.0484570
0	Inf	1	r.squared	0.9999997

```
summary(my.results.automatic)
```

start	end	auclast	cmax	tmax	half.life	aucinf
0	24	74.6 [24.3]	.	.	.	.
0	Inf	.	8.65 [17.0]	1.14 [0.630, 3.55]	8.18 [2.12]	115 [28.4]

```
my.results.manual <- pk.nca(my.data.manual)
knitr::kable(head(my.results.manual$result))
```

start	end	Subject	PPTESTCD	PPORRES
0	Inf	1	auclast	147.23475
0	Inf	1	cmax	10.50000
0	Inf	1	tmax	1.12000
0	Inf	2	auclast	88.73128
0	Inf	2	cmax	8.33000
0	Inf	2	tmax	1.92000

```
summary(my.results.manual)
```

```
## Warning in summary.PKNCAresults(my.results.manual): No results to summarize
## for aucinf.obs in result row 1
```

start	end	auclast	cmax	tmax	aucinf.obs
0	Inf	98.7 [22.5]	8.65 [17.0]	1.14 [0.630, 3.55]	.